


From: Reagan, Chuck Chuck.Reagan@mksinst.com 
Subject: RE: Ophir 300PDR and Juno
Date: October 4, 2021 at 1:47 PM
To: Warner Miller wam@fau.edu
Cc: Joe Dodd doddj2018@fau.edu

CR

EXTERNAL EMAIL : Exercise caution when responding, opening links, or opening attachments.

Warner,

Thank you very much for the reply and I would love visiting again. Do I need to quote you anything new, or do you already own what you need?

Chuck Reagan
Southeastern Manager
Ophir-Spiricon Products



MKS Instruments, Inc.
3050 North 300 West
North Logan, Utah 84341
Cell: 813-716-8972
Main: 435-753-3729
Sales: 435-755-5421
Fax: 435-753-5231

Chuck.Reagan@mksinst.com

www.ophiropt.com/photonics

www.ophiropt.com

<http://www.mksinst.com>



NEW! ARIEL - SMALL SIZE MEETS HIGH POWER

**AN INDUSTRIAL ALL-IN-ONE SENSOR FOR HIGH POWER
LASER MEASUREMENT**

[LEARN MORE](#)



From: Warner Miller <wam@fau.edu>
Sent: Monday, October 4, 2021 12:10 PM
To: Reagan, Chuck <Chuck.Reagan@mksinst.com>
Cc: Joe Dodd <doddj2018@fau.edu>
Subject: [EXTERNAL] Re: Ophir 300PDR and Juno

This email originated outside of MKS. Use caution when sharing information or opening attachments and links.

Hi Chuck,
Perfect. I will read this over and place my order. We should get you back to FAU for another visit sometime (Winter is a great time!). I am Cc'ing my student (Joe Dodd).
Thank you

Thank you.
Best Wishes,
Warner

Dr. Warner A. Miller
Professor of Physics
Department of Physics
Florida Atlantic University
wam@fau.edu
<http://www.physics.fau.edu/~wam>

—
{0,0}
l)___)

—''—''—

On Oct 4, 2021, at 12:52 PM, Reagan, Chuck
<Chuck.Reagan@mksinst.com> wrote:

EXTERNAL EMAIL : Exercise caution when responding, opening links, or opening attachments.

Warner,

One of my guys answered me with the following notes. See also the attached files.

Linux and Raspberry PI **Linux and Raspberry PI**

Linux Project:

Attached is the Linux Eclipse project you can also download the package from:

<https://drive.google.com/file/d/1gcIJBQ-tlhrrva64kSCoYOP22wKmVcZz/view?usp=sharing>

This package enables customers to integrate communication with Ophir [USB talking meters](#), or [USB PC interfaces](#) such as the Juno+, via our command set – see the attached pdf document.

In the main.cpp file there is sample code of how to work with the library.

There are two ways of receiving the measurement value:

There are two ways of receiving the measurement value.

- Command mode: e.g. "SP"- Send Power command (polling mode) via `usbDriver.executeCommand`
- Continuous mode (streaming measurement mode): via `usbDriver.startMeasuring` and then retrieve the measurement by calling to `usbDriver.getReading`

The following Ophir instruments are presently supported in `Definition.h`:

```
enum SupportedProducts {  
    Nova2 = 0x333,  
    Vega = 0x334,  
    StarLite = 0x345,  
    StarBright = 0x346,  
    Juno = 0x777,  
};
```

Compatibility with other meters and relatively newer meters and interfaces, the Linux code needs to be modified to add the PID's of the required instrument, and the streaming commands. The following tables indicates PIDs and which existing instrument to copy streaming commands from.

Instrument	PID	Streaming commands like
USBI	0222	Novall/Vega
Novall	0333	Novall/Vega
Vega	0334	Novall/Vega
StarLite	0345	StarBright
843-R	E345	StarBright
StarBright	0346	StarBright
1919-R	E346	StarBright
PM841PE	E501	Juno
PM844PE	E502	Juno
Juno	0777	Juno
Juno_Plus	0502	Juno
Centauri	0790	StarBright
x938	E794	StarBright

To support streaming mode, using `usbDriver.startMeasuring`, there is the following list in `UsbDriver.cpp`:

```
mSupportedProductsSet[Nova2] = CsCommand("CS 1 1 3", "CS 0");  
mSupportedProductsSet[Vega] = CsCommand("CS 1 1 3", "CS 0");  
mSupportedProductsSet[Juno] = CsCommand("CS 3", "CS 1");  
mSupportedProductsSet[StarLite] = CsCommand("CS 2", "CS 1");
```

```
mSupportedProductsSet[StarBright] = CsCommand("CS 2", "CS 1");
```

This also needs to be extended to support new instruments. The definitions should be copied from the existing definitions according to the table above.

Here are the new lines. You can amend for just a single instrument such as the Juno+, just make sure to do the same thing in both files.

In Definition.h:

```
enum SupportedProducts {  
    Nova2 = 0x333,  
    Vega = 0x334,  
    StarLite = 0x345,  
    StarBright = 0x346,  
    Juno = 0x777,  
    Pm843 = 0xE345,  
    Pm1919 = 0xE346,  
    JunoPlus = 0x502,  
    Pm841 = 0xE501,  
    Pm844 = 0xE502,  
    Centauri = 0x790,  
    X938 = 0xE794,  
};
```

In UsbDriver.cpp:

```
mSupportedProductsSet[Nova2] = CsCommand("CS 1 1 3", "CS 0");  
mSupportedProductsSet[Vega] = CsCommand("CS 1 1 3", "CS 0");  
mSupportedProductsSet[Juno] = CsCommand("CS 3", "CS 1");  
mSupportedProductsSet[StarLite] = CsCommand("CS 2", "CS 1");  
mSupportedProductsSet[Pm843] = CsCommand("CS 2", "CS 1");  
mSupportedProductsSet[Pm1919] = CsCommand("CS 2", "CS 1");  
mSupportedProductsSet[JunoPlus] = CsCommand("CS 3", "CS 1");  
mSupportedProductsSet[Pm841] = CsCommand("CS 3", "CS 1");  
mSupportedProductsSet[Pm844] = CsCommand("CS 3", "CS 1");  
mSupportedProductsSet[Centauri] = CsCommand("CS 2", "CS 1");  
mSupportedProductsSet[X938] = CsCommand("CS 2", "CS 1");
```

Raspberry Pi Info:

The following info was provided by a customer noting changes to the Linux project to work with Raspberry PI. Several customers have successfully used these instructions:

The basic steps

the basic steps:

- `sudo apt install libusb-1.0-0-dev`
- `rm include directory`
- modify everywhere `#include <libusb.h>` to `#include <libusb-1.0/libusb>`
- added `-lpthread` to line 7 of `Default/objects.mk`
- execute `make` in `Default` directory
- execute binary as superuser (`sudo`)

Linux Project:

Attached is the Linux Eclipse project you can also download the package from:

<https://drive.google.com/file/d/1gclJBQ-tlhrrva64kSCoYOP22wKmVcZz/view?usp=sharing>

This package enables customers to integrate communication with Ophir USB talking meters, or USB PC interfaces such as the Juno+, via our command set – see the attached pdf document.

In the `main.cpp` file there is sample code of how to work with the library.

There are two ways of receiving the measurement value:

- Command mode: e.g. “SP”- Send Power command (polling mode) via `usbDriver.executeCommand`
- Continuous mode (streaming measurement mode): via `usbDriver.startMeasuring` and then retrieve the measurement by calling to `usbDriver.getReading`

The following Ophir instruments are presently supported in `Definition.h`:

```
enum SupportedProducts {  
    Nova2 = 0x333,  
    Vega = 0x334,  
    StarLite = 0x345,  
    StarBright = 0x346,  
    Juno = 0x777,  
};
```

Compatibility with other meters and relatively newer meters and interfaces, the Linux code needs to be modified to add the PID's of the required instrument, and the streaming commands. The following tables indicates PIDs and which existing instrument to copy streaming commands from.

Instrument	PID	Streaming commands like
USBI 0222	Novall/Vega	
Novall 0333	Novall/Vega	
Vega 0334	Novall/Vega	
StarLite 0345	StarBright	
843-R E345	StarBright	
StarBright	0346	StarBright
1919-R E346	StarBright	
PM841PE	E501	Juno
PM844PE	E502	Juno
Juno 0777	Juno	
Juno_Plus	0502	Juno
Centauri	0790	StarBright
x938 E794	StarBright	

```
mSupportedProductsSet[Nova2] = CsCommand("CS 1 1 3", "CS 0");  
mSupportedProductsSet[Vega] = CsCommand("CS 1 1 3", "CS 0");  
mSupportedProductsSet[Juno] = CsCommand("CS 3", "CS 1");  
mSupportedProductsSet[StarLite] = CsCommand("CS 2", "CS 1");  
mSupportedProductsSet[Pm843] = CsCommand("CS 2", "CS 1");  
mSupportedProductsSet[Pm1919] = CsCommand("CS 2", "CS 1");  
mSupportedProductsSet[JunoPlus] = CsCommand("CS 3", "CS 1");  
mSupportedProductsSet[Pm841] = CsCommand("CS 3", "CS 1");  
mSupportedProductsSet[Pm844] = CsCommand("CS 3", "CS 1");  
mSupportedProductsSet[Centauri] = CsCommand("CS 2", "CS 1");  
mSupportedProductsSet[X938] = CsCommand("CS 2", "CS 1");
```

Raspberry Pi Info:

The following info was provided by a customer noting changes to the Linux project to work with Raspberry PI. Several customers have successfully used these instructions:

The basic steps

- sudo apt install libusb-1.0-0-dev
- rm include directory
- modify everywhere #include <libusb.h> to #include <libusb-1.0/libusb>
- added -lpthread to line 7 of Default/objects.mk
- execute make in Default directory
- execute binary as superuser (sudo)

Chuck Reagan
Southeastern Manager
Ophir-Spiricon Products
<image001.png>
MKS Instruments, Inc.
3050 North 300 West
North Logan, Utah 84341
Cell: 813-716-8972
Main: 435-753-3729
Sales: 435-755-5421
Fax: 435-753-5231

Chuck.Reagan@mksinst.com
www.ophiropt.com/photonics

www.ophiropt.com
<http://www.mksinst.com>
<image007.jpg>

<image004.png>

From: Reagan, Chuck
Sent: Monday, October 4, 2021 9:49 AM
To: wam@fau.edu
Subject: Ophir 300PDR and Juno

Warner,

Thank you for your interest in our products. I have sent your request to my power meter guru. He will have to consult with the factory engineers to see if this will work. Give me a few days to collect the answers. It has been many years since I last visited your campus.

Name:	Warner Miller
Company:	Florida Atlantic University
State:	Florida
Country:	United States
E-Mail:	wam@fau.edu
Message:	I am looking to put a compact power meter coupling into either an SMA or FC fiber for low power (laser pointer). It has to be light weight as I will be putting this on a gimbal. I would also like to interface this to a raspberry pi 4 computer running linux. I was thinking of the 300PDR connected to the Juno. What do you recommend? Can this run under linux? Thank you.
Page:	PD300R

Chuck Reagan
Southeastern Manager
Ophir-Spiricon Products
<image001.png>
MKS Instruments, Inc.
3050 North 300 West
North Logan, Utah 84341
Cell: 813-716-8972
Main: 435-752-2720

main: 435-753-5729
Sales: 435-755-5421
Fax: 435-753-5231

Chuck.Reagan@mksinst.com

www.ophiropt.com/photonics

www.ophiropt.com

<http://www.mksinst.com>

<image008.jpg>

[<image004.png>](#)

<ophiropt.drivers.usb.linux.zip><Ophir User Commands.pdf>